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EN INSTRUCTION MANUAL

1 Specific purpose

BetaBASE® is a bioresorbable bone replacement material for solid tissue replacement. It serves for the filling of mechanically largely non-stressed, pathological, traumatic or surgery-related bone defects in orthopaedics/traumatology.

Before using it is recommended to verify the suitability of the products for the planned intervention.

Indications

- As addition to autogenic or allogenic spongiosa with spondylodesis, with vertebral body replacement and with reconstructions of bone defects
- Filling of defects with correction osteotomies
- Filling of defects in bone traumatology
- Filling of defects after tumour resection
- Application with arthrodeses
- Filling of the extraction locations of autogenic spongiosa

Contraindications

- Medication and therapies that influence bone wound healing, e.g.: Diabetes mellitus, treatments with glucocorticoids and antineoplasm acting agents or irradiation with γ -radiation
- Osteomyelitis
- Infections at the site of the implantation

2. Factors that can impair the success of the operation

- Insufficient bone quality (e.g. osteoporosis)
- Local bone tumours
- Systemic diseases and metabolic disorders
- Anamnestic indications for infections and falls
- Drug dependency, alcohol and medication abuse
- Mental inability of the patient to understand and follow the instructions given by the physician

Possible negative effects

With the application of BetaBASE® it can in rare cases come to the occurrence of the following product, respectively procedure related negative effects.

- Infections
- Bone fractures or pseudarthrosis of an existing fracture
- Bleeding events
- Delayed/disturbed wound healing
- Injuries that are caused through the implantation, including vascular injuries, nerve lesions and tissue injuries
- Heterotopic ossifications
- Implant material intolerance reactions

3. Patient briefing

The physician is to inform the patient about possible impairments of the operation success and possible negative effects according to the sections “contraindications” and “possible negative effects”, insofar as applicable for the patient. Furthermore, the patient is to be informed, which measures he or she can undertake to reduce possible impairments and negative effects.

4. Operation technique

A surgeon that is not familiar with the specific product is expected to exactly study the surgery technique before the operation.

5. Handling

NOTE: This product is packed in a protective and sterile packaging.
WARNING: Danger of injury and risk of damage of products: The non-observance of these instructions for use can lead to injuries to patients, users and third parties, as well as to damage of the product. Carefully read through the instructions for use and always adhere to the described instructions. Verify the function of the products used in combination.
WARNING: Before using verify the expiry date, as well as the packing for any damage. Products with overdue expiry date, respectively damaged packaging must not be used.
WARNING: Risk of infection: Through the use of non-sterile implants there is risk of infection for the patient. Only remove the implants from the sterile packing immediately before the application.
WARNING: Risk of infection: A reconditioning of the disposable article can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

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EN INSTRUCTION MANUAL

7. Application

BetaBASE® should exclusively be used by a skilled user with corresponding training and experience in orthopaedics/traumatology. It is recommended to use the granulate in moist condition, the mass of the individual particles adheres in moist condition and allows more accurate insertion into the defect than in dry form. For moistening, sterile isotonic saline solution can be used, as well as own blood of the patient or blood plasma (PRP), with suspected infections under addition of an antibiotic solution. If the defect is larger than 2 cm³ we recommend mixing BetaBASE® with autologous bone material, if the removal is anatomically possible. During the application, the material should not be compressed but evenly layered into the defect, respectively filled to maintain the interconnecting porosity. The selection of the particle size to be used generally depends on the size of the defect. Thereby the following rule of thumb applies: The larger the size of the defect, the larger the selected particle size should be. Depending on the size, the introduction into the defect can be performed with a spoon or spatula. After the filling of the bone defect, which should take place volume filling up to the edges, the wound is closed with the usual suturing technique. Depending on the regeneration dynamics of the patient, the resorption time is approximately 6 to 24 months.

8. Traceability

The self-adhesive identification labels that are enclosed with the original packing must be stored for traceability and be adhered to the operation report, as well as in the patient file. The individual responsible for the stock-keeping must maintain records of the available products, including article number (REF) and batch number (LOT), so that in the event of changes or recalls action is quickly possible.

9. Packaging and sterility

Every individual product is double-packed and sterilised through gamma radiation with a dose of at least 25 kGy. The sterility is guaranteed up to the expiry date, insofar as the packing is neither opened nor damaged.

10. Product material

BetaBASE® is a pure β -tricalcium phosphate ceramic. In the individual granulate particles there are micro pores (up to 5 μ m), as well as also macro pores (up to 1 mm).

Composition:
Calcium oxide CaO 52.0 to 54.2 weight %
Phosphorus pentoxide P₂O₅ 45.8 to 48.0 weight %
The ratio Ca/P is 1.5.

11. Reusability

WARNING: Risk of infection: A reconditioning of the disposable article can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

12. Storage conditions

The product must be stored dry in the unopened original packing. Products with elapsed expiry date, respectively damaged packing must not be used.

13. Responsibility

Complications or other negative effects due to incorrect indication or surgical technique, unsuitable selection or handling of products or instruments, asepsis, etc. fall under the responsibility of the surgeon.

14 Conformity with directives

According to the medical device directive (MDD) 93/42/EEC this medical device is labelled with the CE label. If there is an identification number after the CE label then this designates the responsible specified authority.

15 Disposal

No special measures are necessary during the disposal of this medical device. The national laws and directives must be observed.

16 Packing symbols

	Manufacturer
	Observe instructions for use
	Use before
	Do not use product with damaged packing
	Sterilised through irradiation

	Do not reuse
	Do not re-sterilise
	Batch number
	Article number
	Store dry

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