
EpiGARD®

CE 0297

BIOVISION

BIOMATERIAL

BIOVISION GmbH
Am Vogelherd 52
98693 Ilmenau
Germany

info@biovision.de
+49 3677 64 07 0
+49 3677 64 07 13

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EN INSTRUCTION MANUAL

1. Intended use

EpiGARD® is a synthetic skin replacement for temporary wound covering, for wound cleaning and the conditioning of superficial wounds.

Indications

With external wounds with open wound treatment and for preparation of a secondary wound closure, such as:

- Superficial defect wounds
- Wounds after surgical necrosis removal with burns
- Operation wounds
- Wounds at removal sites for skin transplants

Contraindications

- Wounds with pocket formation
- Wound infections
- Severe wound secretion
- Bleeding tendency
- Simultaneous application of creams or ointments

2. Possible negative effects

With exchange of the wound covering patient-dependent pain and slight bleeding can occur. The specialist user should undertake individual measures for pain reduction. A renewal of the wound covering after a short interval can also significantly reduce the pain.

There is a potential risk of the wound becoming infected during application or when changing the dressing.

3. Patient briefing

The specialist user is to inform the patient about possible negative effects according to the sections "Contraindications" and "Possible negative effects", insofar as applicable for the patient. Furthermore, the patient is to be informed, which measures he or she can undertake to reduce possible impairments and negative effects.

4. Handling

NOTE: This product is packed in a protective and sterile packaging.

WARNING: Danger of injury: The non-observance of these instructions for use can lead to injuries of patients, users and third parties, as well as to damage of the product. Carefully read through the instructions for use and always adhere to the described instructions. Verify the function of the products possibly used in combination (suture material, dressing material, ...).

WARNING: Before using verify the expiry date as well as the packing for any damages. Products with elapsed expiry date, respectively damaged packing must not be used.

WARNING: Risk of infection: Reconditioning of the disposable article can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

Re-sterilisation is not permitted.

5. Variants

EpiGARD® is available as follows:

Article number	Size	Contents of the packaging
WV3200	8 x 5 cm	10 Pcs/package unit
WV3201	8 x 10 cm	10 Pcs/package unit
WV3202	8 x 23 cm	10 Pcs/package unit
WV3203	12 x 30 cm	10 Pcs/package unit
WV3204	12 x 45 cm	10 Pcs/package unit

6. Application

EpiGARD® may only be applied through correspondingly qualified specialist personnel. Before using EpiGARD® the wound area must be carefully cleaned. Ointment residues and larger necroses are to be removed.

Place EpiGARD® with the open pore foam layer onto the wound!

To ensure a close contact with the wound surface, EpiGARD® must be trimmed according to the wound size, under observation of sterile conditions and not extend beyond the wound edges. It should be held in position e.g. with a light bandage. Due to frequent change of dressing saturation should be dispensed with.

EpiGARD® is to be regularly examined and with fluid accumulation or irritation of the wound edges it is to be replaced.

In general daily renewal of EpiGARD® is recommended.

After every exchange a decision should be undertaken about the definitive wound closure. EpiGARD must not be applied cumulatively for longer than 30 days. With contaminated wounds such as e.g. with ulcerations, residual necroses and wound exudates can be removed from the wound area through frequent change of dressing (1 to 2 times daily). The removal of the infection is the requirement for the adhesiveness of EpiGARD® on the wound and the granulation stimulation.

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EN INSTRUCTION MANUAL

7. Packaging and sterility

Each individual product is double packed and sterilised by means of ethylene oxide. The sterility is guaranteed up to the expiry date, insofar as the packing is neither opened nor damaged.

8. Product material

EpiGARD® is a nontextile, micro porous plastic material without medication impregnation. The material is double layered. The lower side forms an open matrix made of elastic soft polyurethane foam. The other side consists of a thin, micro porous polytetrafluorethylene film, is permeable to air and water vapour, however prevents the migration of bacteria and the passage of fluid and secretion.

9. Reusability

WARNING: Risk of infection: A reconditioning of the disposable article can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

10. Storage conditions

The product must be stored dry in the unopened original packing. Products with elapsed expiry date, respectively damaged packing must not be used.

11. Responsibility

Complications or other negative effects due to incorrect indication, unsuitable selection or handling of products, asepsis etc. fall under the responsibility of the user.

12. Conformity with directives

According to the medical device directive (MDD) 93/42/EEC this medical device is labelled with the CE label. If there is an identification number after the CE label then this designates the responsible specified authority.

13. Disposal

EpiGARD® must be disposed of appropriately. The national laws and directives must be observed.

14. Packing symbols

	Manufacturer
	Observe instructions for use
	Use by
	Do not use with damaged packing
	Sterilised with ethylene oxide
	Do not reuse
	Do not re-sterilise
	Store dry
	Protect from sunlight

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