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## 1. Intended use

The PolyPIN® is a bioresorbable pin for bone or osteochondral fragment fixation. Before using it is recommended to check the suitability of the products for the planned intervention.

### Indications

#### PolyPIN® 1.5

- Bone or osteochondral fragment fixation, e.g. lateral clavicle fractures, metacarpal head fractures, fragment attachment
- Arthrodeses of fingers and toes
- Stabilisation of certain finger fractures
- Osteochondral fractures and dissecans

#### PolyPIN® 2.0

- Apical fragments, e.g. fractures of the radius head, patella edge fractures, fractures at the proximal and distal end of the Ossa metatarsalia and metacarpalia
- Osteochondral fractures or dissecans, e.g. talar dome and femur condyle
- Spongiosa or lightly load-bearing fragments, e.g. layer-wise restructuring of calcaneus or acetabulum fractures, cortico spongiosa fragment transplants

#### PolyPIN® 2.7

- Large osteochondral fragments, e.g. hip head fractures (Pipkin fractures), malleolar fractures (Weber A fractures), radius fractures, longitudinal patella fractures

### Contraindications (for all sizes of PolyPIN®)

- Patients with pronounced osteoporosis
- Infantile fractures in the area of the epiphyseal cartilage
- Fractures with high biomechanical load-bearing (e.g. load-bearing osteosyntheses with supracondylar humerus fractures)
- Medications and therapies influencing bone wound healing
- Osteomyelitis
- Infections at the implantation site

## 2. Factors that can impair the success of the operation

- Insufficient bone quality (e.g. osteoporosis)
- Local bone tumours
- Systemic diseases and metabolic disorders
- Anamnestic indications for infections and falls
- Drug dependency, alcohol and medication abuse
- Mental inability of the patient to understand and follow the instructions given by the physician
- Pin head dislocations or damage to a joint surface opposite to the fracture location through an incomplete counter sinking of the PolyPIN® head
- Fragment dislocation through incorrect indication selection or too early, uncontrolled post operative load exertion
- A continuous ligament and tendon traction can lead to failed interventions through secondary dislocation

### Possible negative effects

In rare cases, the following product- or procedure-related negative impacts can occur with the use of PolyPIN®:

- Breakage of the PolyPIN®
- Foreign body reactions like local edemata and osteolysis
- Infections
- Bleedings
- Limited mobility
- Heterotopic ossifications
- Complex regional pain syndrome (CRPS)
- Malpositions like Cubitus varus
- Dislocation and migration
- Nerve injury

## 3. Patient information briefing

The physician is to inform the patient about possible impairments of the operation success and possible negative effects according to the sections "contraindications" and "possible negative effects", insofar as they are applicable for the patient. Furthermore, the patient is to be informed, which measures he or she can undertake to reduce possible impairments and effects.

The patient must be informed that the product materials can cause allergic reactions including reactions to foreign bodies, tissue irritation and inflammation. If sensitivity to foreign bodies is suspected, a suitable test should be performed prior to implantation. If the implant has to be removed, the potential risk of a second surgical procedure must be taken into account. Appropriate post-operative care should be administered after the implant is removed.

## 4. Operation technique

For the successful use of the PolyPIN®, it is necessary to properly apply pre-operative and operative procedures and surgical methods, and to ensure that the appropriate implant is selected and suitably positioned. The surgeon must select the suitable size of article in accordance with the diagnosis and surgical method used.

**WARNING:** The surgeon must explain appropriate post-care and rehabilitation measures to the patient so the implant is not prematurely exposed to loading and to prevent other complications.

A: Locking screw, B: Retaining screw, C: Shearing head, D: Limit stop, E: PolyPIN® head, F: Sharpener

### ① Making the borehole

After the fragment is repositioned, a borehole must be inserted as perpendicular as possible to the surface of the fracture, with a diameter equal to that of the selected PolyPIN®.

PolyPIN® 1.5	1.5 mm borehole
PolyPIN® 2.0	2.0 mm borehole
PolyPIN® 2.7	2.7 mm borehole

If more than one PolyPIN® is used for the fixation procedure, the boreholes should converge rather than be parallel to one another. This ensures that the fragment is securely positioned.

### ② Measuring length

After the locking screw (A) has been loosened, the measuring point of the length measuring device must be inserted down to the bottom of the borehole and the depth of the bore determined. Once this has been performed, the locking screw (A) is tightened. This will transfer the measured length directly to the cutting gauge fixed in place to the rear. After the subsequent shearing, the PolyPIN® will be approx. 1 mm shorter than the measured borehole depth.

### ③ Shortening the PolyPIN®

The PolyPIN® is firstly inserted head first into the measuring device until it reaches the limit stop at the end. It must then be secured in place with the retaining screw (B), with the shearing head placed on it so that the arrow sign is visible. Turning the shearing head (C) in the direction of the arrow will cut the PolyPIN® to the desired length.

Warning: Do not turn the shearing head against the direction of the arrow, as this can press the PolyPIN® against the retaining screw (B) and cause it to become deformed

### ④ Sharpening the point of the PolyPIN®

For easier insertion, the pin is sharpened with the sharpener (F). To do this, after removing the shearing element, the pin should be ejected from the length measuring device and then securely fixed in place again using the retaining screw (B). Then use the sharpener (F) to gently sharpen the PolyPIN® in the same way as sharpening a pencil.

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## 8. Application

**ATTENTION:** The product must be exclusively implanted with the instruments intended for this purpose. During the implantation the patient must be treated and observed with the usual medical diligence. This includes the control of the progression of the procedure, and the monitoring of the vital parameters and the anaesthesia.

## 9. Traceability

The self-adhesive identification labels that are enclosed with the original packing must be stored for traceability and be adhered to the operation report, as well as the patient file. The individual at the hospital responsible for the stock keeping must maintain records of the available products, including article number (REF) and batch number (LOT), so that in the event of changes or recalls action can be quickly taken.

## 10. Packaging and sterility

Every individual product is double packed and sterilised through gamma radiation with a dose of at least 25 kGy. The sterility is guaranteed up to the expiry date, insofar as the packing is neither opened nor damaged.

## 11. Product material

The bioresorbable PolyPIN® is made of poly(L-co-D,L-lactate 70/30). This bioresorbable material is biocompatible and degradable.

## 12. Reusability

**WARNING:** Risk of infection: Reconditioning of disposable articles can endanger the functionality and safety of the product and lead to infections. You must not recondition disposable articles. Dispose of disposable articles appropriately.

## 13. Storage conditions

The product must be stored dry at max. 25 °C in the unopened original packing. Products with elapsed expiry date, respectively damaged packing must not be used.

## 14. Temperature indicator

PolyPIN® implants may not be exposed to temperatures in excess of 40 °C (104 °Fahrenheit), even for a brief time. The implant can become deformed if exposed to excessive heat. If the product was ever exposed to excessive heat, the colour of the indicator incorporated into the outer packaging will change from light grey to black. It is prohibited to use products with an indicator that is coloured black!

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## 15. Responsibility

Complications or other negative effects due to incorrect indication or surgical technique, unsuitable selection or handling of products or instruments, asepsis etc. fall under the responsibility of the surgeon.

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## 16. Conformity with directives

This medical device is marked with CE labelling according to the Medical Device Directive (MDD) 93/42/EEC. If there is an identification number after the CE labelling then this designates the responsible specified authority.

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## 17. Disposal

No special measures are necessary for the disposal of this medical device. The national laws and directives must be observed.

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## 18. Packing symbols

	Manufacturer		Do not use the product with damaged packing		Do not resterilise
	Observe the instructions for use		Sterilised through irradiation		Temperature limit
	Maximum shelf life of the product		The product is not suitable for reuse		Store dry

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